

EVALUATION OF HCV RNA TESTING FROM FINGER-STICK CAPILLARY DRIED BLOOD SPOT AND VENEPUNCTURE-COLLECTED SAMPLES: A COHORT STUDY

Catlett B^{1,2}, Bajis S², Starr M¹, Dore GJ², Cunningham PH^{1,2}, Hajarizadeh B², Grebely J² and Applegate TL^{2*}

¹ NSW State Reference Laboratory for HIV, St Vincent's Centre for Applied Medical Research, Darlinghurst, Sydney, Australia

² The Kirby Institute, UNSW Sydney, Sydney, New South Wales, Australia

*Joint senior authors, equally contributed to the work

Background: Simplified HCV diagnostic strategies increasing testing and linkage to care are needed. Novel collection methods e.g. finger-stick capillary dried blood spots (DBS) have advantages over standard phlebotomy, especially among people who inject drugs where venous access may be a barrier to testing. The study aim was to evaluate diagnostic performance of the Hologic Aptima HCV Quant assay for detection and quantification of HCV RNA with finger-stick capillary DBS and venepuncture samples.

Methods: Plasma and finger-stick capillary DBS samples were collected from participants in an observational cohort enrolled at six sites in Australia. We compared sensitivity, specificity, bias and agreement of the Aptima HCV assay for HCV RNA detection by finger-stick capillary DBS with venepuncture (gold standard). Retesting was performed on samples that had undetectable HCV RNA in plasma and detectable/unquantifiable HCV RNA (<10 IU/mL) in DBS.

Results: Of 175 participants enrolled in December 2016, 164 participants (11 participants on DAA therapy excluded) had paired viral load results. HCV RNA was detected in 45 (27%) of 164 participants based on the gold standard. Among 38 people with undetectable HCV RNA in plasma and detectable/unquantifiable HCV RNA (<10 IU/mL) in DBS, 32 samples were available for re-testing (all HCV undetectable). Sensitivity for HCV RNA detection in DBS was 95.6% (95% CI 84.9-99.5%), specificity was 94.1% (95% CI 88.3-97.6%). Sensitivity for HCV RNA quantification in DBS (≥ 10 IU/mL in plasma) was 100% (95% CI 91.8%-100%), specificity was 100% (95% CI 97%-100%). Finally, a small bias of 0.37 Log₁₀ in plasma over DBS (95% CI 0.78-1.55) was observed with good agreement ($R^2=0.96$).

Conclusion: The Aptima HCV assay can detect active infection from DBS with acceptable diagnostic performance and is clinically comparable to plasma. Repeat testing of HCV RNA detectable/unquantifiable samples suggest that a quantifiable result is a more reliable indicator for active infection.

Disclosure of Interest Statement:

The views expressed in this publication do not necessarily represent the position of the Australian government.

Hologic provided Aptima HCV Dx reagents. The Kirby Institute is funded by the Australian Government Department of Health and Ageing and the National Health and Medical Research Program Grant (#1053206). B.C and SB are supported by an Australian Postgraduate Award from the University of New South Wales (UNSW) Sydney and an Australian Government Research Training Program Scholarship. J.G

is supported by a National Health and Medical Research Council Career Development Fellowship.

T.A has received travel grants and travel support from Abbott Diagnostics and travel grants from Cepheid. J.G is a consultant/advisor and has received research grants from AbbVie, Bristol-Myers Squibb, Cepheid, Gilead Sciences and Merck/MSD.