


ViziShot 2

For Improved Access and Enhanced
Lymph Node Penetration



ViziShot 2


Improved Access and Enhanced Lymph Node Penetration for More Efficient Sampling in EBUS-TBNA

Efficient Echogenicity
for clear ultrasound visibility.



16153

Green Color for Increased Sheath Visibility
supporting visualization during endoscopic observations.



Ergonomic Handle Design
supporting better needle control during the procedure.




Rotatable Sheath Adjuster
for convenient adjustment of the sheath length.



Improved Handle Fixation
allowing firm locking on the scope.





Sharper Needle Tip Design
enhancing puncture capability and supporting smooth penetration of the bronchial wall.



Improved Flexibility
enabling easier access to challenging lymph nodes and stations.

ViziShot 2

Specifications

EBUS-TBNA Needles

ViziShot 2

Article Name	Min. Working Channel Ø	Needle Gauge	Max. Needle Length	Quantity
NA-U401SX-4021	2.0 mm	21G	40 mm	5 sets
NA-U401SX-4022	2.0 mm	22G	40 mm	5 sets

ViziShot 2 FLEX

Article Name	Min. Working Channel Ø	Needle Gauge	Max. Needle Length	Quantity
NA-U403SX-4019	2.2 mm	19G	40 mm	1 set

ViziShot 2 25G

Article Name	Min. Working Channel Ø	Needle Gauge	Max. Needle Length	Quantity
NA-U401SX-4025N	2.0 mm	25G	40 mm	5 sets

ViziShot

Article Name	Min. Working Channel Ø	Needle Gauge	Max. Needle Length	Quantity
NA-201SX-4021	2.0 mm	21G	40 mm	5 sets
NA-201SX-4022	2.0 mm	22G	40 mm	5 sets

www.olympus.eu/et-catalog



As medical knowledge is constantly growing, technical modifications or changes of the product design, product specifications, accessories and service offerings may be required.

OLYMPUS

OLYMPUS EUROPA SE & CO. KG
Postbox 10 49 08
20034 Hamburg, Germany
www.olympus-europa.com