




# BONASTENT®

BONA Opens New Area!



Standard Sci-Tech Inc.


Polygon Mesh Surface Structure




Adapting Flexibility




Segmental Compression



Low Returning Force



Repositionable



## The Difference?

A Company



Hook Wire Structure

- Difficult Penetration due to Larger Profile

B Company



Cross Wire Structure

- Inaccurate Positioning due to Greater Shortening & Migration

BONASTENT®



Hook & Cross Wire Structure (Polygon Mesh Surface)

- Low Rate of Shortening for Accurate Positioning
- Small Cavities for In-Growth Reduction
- Low Returning Force for Conformability
- Segmental Compression for Migration Reduction
- Reduced Delivery Device Diameter

# BONA Opens New Area!

## Product Catalog

S t a n d a r d S c i . T e c h I n c .

Esophageal & Esophageal with SHIM Anti-Reflux Valve  
Biliary  
M-Hilar  
M-Intraductal  
AL, Hybrid & DEUS  
Duodenal / Pyloric  
M-Duodenal  
Colo-Rectal  
Tracheal / Bronchial

## BONASTENT®

Since our first product was introduced in the market in 2005, we have been promoting our products globally and positioned ourselves as one of the top stent producers.

Hurdling over the barriers, we have been growing steadily in the competitive international market with high quality features and test proven performance results, and the global users, based on their experience, have gladly accepted our proposal to add BONASTENT® on their supply list.

We fully acknowledge that our success was only possible with your complete support and precious comments for us to strive for further improvements, and we will continue to bear in mind that our business is closely connected with all of you.

Look forward to our new products in the market, the results of our endless innovation and creativity, and we sincerely hope to make that connection with you soon.

BONASTENT® Tracheal/Bronchial

# Tracheal/Bronchial




## Features

- Silicone covering
- Non-irritating ends
- Segmental compression
- Repositionable
- Platinum radiopaque markers

# Ordering Information

## BONASTENT® Tracheal/Bronchial



| Catalog No. | Stent         |                   | Delivery Device |                    |
|-------------|---------------|-------------------|-----------------|--------------------|
|             | Diameter [mm] | Total Length [mm] | Diameter [Fr]   | Usable Length [mm] |
| BTB-060209  | 6             | 20                | 8 (2.66mm)      | 900                |
| BTB-060309  | 6             | 30                | 8 (2.66mm)      | 900                |
| BTB-060409  | 6             | 40                | 8 (2.66mm)      | 900                |
| BTB-060509  | 6             | 50                | 8 (2.66mm)      | 900                |
| BTB-060609  | 6             | 60                | 8 (2.66mm)      | 900                |
| BTB-060709  | 6             | 70                | 8 (2.66mm)      | 900                |
| BTB-060809  | 6             | 80                | 8 (2.66mm)      | 900                |
| BTB-080209  | 8             | 20                | 8 (2.66mm)      | 900                |
| BTB-080309  | 8             | 30                | 8 (2.66mm)      | 900                |
| BTB-080409  | 8             | 40                | 8 (2.66mm)      | 900                |
| BTB-080509  | 8             | 50                | 8 (2.66mm)      | 900                |
| BTB-080609  | 8             | 60                | 8 (2.66mm)      | 900                |
| BTB-080709  | 8             | 70                | 8 (2.66mm)      | 900                |
| BTB-080809  | 8             | 80                | 8 (2.66mm)      | 900                |
| BTB-100209  | 10            | 20                | 8 (2.66mm)      | 900                |
| BTB-100309  | 10            | 30                | 8 (2.66mm)      | 900                |
| BTB-100409  | 10            | 40                | 8 (2.66mm)      | 900                |
| BTB-100509  | 10            | 50                | 8 (2.66mm)      | 900                |
| BTB-100609  | 10            | 60                | 8 (2.66mm)      | 900                |
| BTB-100709  | 10            | 70                | 8 (2.66mm)      | 900                |
| BTB-100809  | 10            | 80                | 8 (2.66mm)      | 900                |
| BTB-120209  | 12            | 20                | 9 (3mm)         | 900                |
| BTB-120309  | 12            | 30                | 9 (3mm)         | 900                |
| BTB-120409  | 12            | 40                | 9 (3mm)         | 900                |
| BTB-120509  | 12            | 50                | 9 (3mm)         | 900                |
| BTB-120609  | 12            | 60                | 9 (3mm)         | 900                |
| BTB-120709  | 12            | 70                | 9 (3mm)         | 900                |
| BTB-120809  | 12            | 80                | 9 (3mm)         | 900                |

\* 500mm delivery device is available.

| Catalog No. | Stent         |                   | Delivery Device |                    |
|-------------|---------------|-------------------|-----------------|--------------------|
|             | Diameter [mm] | Total Length [mm] | Diameter [Fr]   | Usable Length [mm] |
| BTB-140209  | 14            | 20                | 12 (4mm)        | 900                |
| BTB-140309  | 14            | 30                | 12 (4mm)        | 900                |
| BTB-140409  | 14            | 40                | 12 (4mm)        | 900                |
| BTB-140509  | 14            | 50                | 12 (4mm)        | 900                |
| BTB-140609  | 14            | 60                | 12 (4mm)        | 900                |
| BTB-140709  | 14            | 70                | 12 (4mm)        | 900                |
| BTB-140809  | 14            | 80                | 12 (4mm)        | 900                |
| BTB-160209  | 16            | 20                | 12 (4mm)        | 900                |
| BTB-160309  | 16            | 30                | 12 (4mm)        | 900                |
| BTB-160409  | 16            | 40                | 12 (4mm)        | 900                |
| BTB-160509  | 16            | 50                | 12 (4mm)        | 900                |
| BTB-160609  | 16            | 60                | 12 (4mm)        | 900                |
| BTB-160709  | 16            | 70                | 12 (4mm)        | 900                |
| BTB-160809  | 16            | 80                | 12 (4mm)        | 900                |
| BTB-180209  | 18            | 20                | 15.6 (5.2mm)    | 900                |
| BTB-180309  | 18            | 30                | 15.6 (5.2mm)    | 900                |
| BTB-180409  | 18            | 40                | 15.6 (5.2mm)    | 900                |
| BTB-180509  | 18            | 50                | 15.6 (5.2mm)    | 900                |
| BTB-180609  | 18            | 60                | 15.6 (5.2mm)    | 900                |
| BTB-180709  | 18            | 70                | 15.6 (5.2mm)    | 900                |
| BTB-180809  | 18            | 80                | 15.6 (5.2mm)    | 900                |
| BTB-200209  | 20            | 20                | 15.6 (5.2mm)    | 900                |
| BTB-200309  | 20            | 30                | 15.6 (5.2mm)    | 900                |
| BTB-200409  | 20            | 40                | 15.6 (5.2mm)    | 900                |
| BTB-200509  | 20            | 50                | 15.6 (5.2mm)    | 900                |
| BTB-200609  | 20            | 60                | 15.6 (5.2mm)    | 900                |
| BTB-200709  | 20            | 70                | 15.6 (5.2mm)    | 900                |
| BTB-200809  | 20            | 80                | 15.6 (5.2mm)    | 900                |

\* 500mm delivery device is available.

# Delivery Devices

BONASTENT®

- \* Superior Pushability and Flexibility
- \* Kinking Resistance
- \* Design for Easy & Accurate Placement

## Y - Shape Handle with Braided Tubing

10Fr or Less Delivery Device for

- Biliary Non-Covered & Covered
- Duodenal / Pyloric Non-Covered & Covered
- Colo-Rectal Non-Covered & Covered



Radiopaque Band

Radiopaque Marking for Distal End of Delivery Device



Yellow Marking Band

Endoscopic View Marking for Proximal End of Captured Stent



## I - Shape Handle with PTFE Tubing

Larger than 10Fr Delivery Device for

- Esophageal & Esophageal with SHIM Anti-Reflux Valve
- Colo-Rectal Non-Covered & Covered
- Tracheal / Bronchial



Centimeter Marking

Positioning Measurement for Esophageal Stent & Esophageal with SHIM ARV Stent

CE  
0197



Obelis S.A

Boulevard Général Wahis, 53, 1030 Brussels, Belgium Phone: +32.2.732.59.54 Fax: +32.2.732.60.03



Manufactured by :

**SEWOON MEDICAL CO., LTD.**

60 Dorim-gil, Ipjang-myeon, Seobuk-gu, Cheonan-si, Chungcheongnam-do, Korea

Phone: +82-41-584-2903~6, +82-2-922-6555 Fax: +82-41-584-2907

[www.sewoonmedical.co.kr](http://www.sewoonmedical.co.kr)



Distributed by :

**Standard Sci-Tech Inc.**

1Fl, Sewoon Medical Bldg., 46 Wangsan-ro, Dongdaemun-gu, Seoul, Korea

Phone: +82-2-922-3605 Fax: +82-2-922-3606

[www.standard4u.com](http://www.standard4u.com)